



8. Adaptors

Adaptors are used when relatively smaller tubes are inserted into the tube cavities. Most adaptors are made from white polyacetal. Since the specific gravity of adaptors is high, the rotation speed is restricted when adaptors are used.



● Name of adaptors



Part No.	Name	Size (cm)	Material	Q'ty
336697A	2S5 Adaptor	1.3 x 4.9	DE	6
336698A	2S13 Adaptor	1.5 x 9.3	DE	6
336699A	3S13 Adaptor	1.5 x 9.3	DE	6
336700A	4S13 Adaptor	1.5 x 9.3	DE	6
336701A	3S16 Adaptor	1.6 x 9.9	DE	6
336702A	4S16 Adaptor	1.6 x 9.9	DE	6
336703A	5S16 Adaptor	1.6 x 9.9	DE	6
336704A	6.5S16 Adaptor	1.6 x 9.9	DE	6
336682A	1.5A12 Adaptor	1.6 x 7.9	DE	10
336683A	2A12 Adaptor	1.6 x 6.3	DE	10
336684A	3A12 Adaptor	1.6 x 6.3	DE	10
336685A	4A12 Adaptor	1.6 x 6.3	DE	10
336686A	5A12 Adaptor	1.6 x 7.9	DE	10
336687A	6.5A12 Adaptor	1.6 x 7.9	DE	10
336688A	4A40 Adaptor	2.6 x 8.5	DE	8
336689A	5A40 Adaptor	2.6 x 8.5	DE	8
336690A	6.5A40 Adaptor	2.6 x 8.5	DE	8
336691A	12A40 Adaptor	2.6 x 8.5	DE	8
336692A	5A94 Adaptor	3.8 x 8.5	DE	6
336693A	6.5A94 Adaptor	3.8 x 8.5	DE	6
336694A	12A94 Adaptor	3.8 x 8.5	DE	6
336695A	40A94 Adaptor	3.8 x 10.1	DE	6

DE:Delrin (Polyacetal)