


Carefully read the centrifuge instruction manual and the rotor instruction manual before using the rotor.

S100AT5 Angle rotor

Max. speed 100,000 rpm
 Max. R.C.F. 550,000 × g
 Rotor nominal capacity 32.0ml = 8 × 4.0ml
 k factor 19



(1) Summary of tubes

Tube				Space cap, Spacer or Crown		Maximum speed (rpm)	Maximum R.C.F. (× g)
Part No.	Name	Size (φ × L cm)	Net capacity (ml)	Part No.	Name		
S404332A	4 PC Tube	1.3×6.0	3.4	—	—	100,000	550,000
S303351A	4 PA Thick-walled tube	1.3×6.0	3.4	—	—	80,000	352,000
345319A	5 PA Seal tube*	1.3×5.3	5.1	S404331	B4-Space cap	100,000	550,000
S308892A	Micro tube 1.5ml	1.09×4.0	1.3	—	—	50,000	107,000
84710901	Sampling tube 1.5ml	1.08×4.2	1.3	—	—	50,000	107,000
(3400-1332)	(3.5ml Ultra-Lok™ tube)	1.3×3.2	3.0	(3405-1003)	(Spacer)	100,000	550,000
(3400-1351)	(5ml Ultra-Lok™ tube)	1.3×5.1	4.1	(3405-1003)	(Spacer)	100,000	550,000
—	1.5ml Micro tube	1.1×4.2	1.0	—	—	50,000	107,000
—	1.5ml Screw cap tube**	1.0×4.6	1.3	—	—	100,000	463,000
S302893A	5 PA Re-seal tube***	1.3×5.1	5.1	S405443A	Crown D2 ass'y	100,000	550,000

Products in parentheses are those of the supply discontinuance.

* Not sold in Europe and U.S.A.

** Applicable tubes

Maker/Brand	Product No.
CORNING	430909
Simport	T334-5/5S

***For Europe and U.S.A.

